

The Soul Code Method Level 2 Procedure Sheet

Client name: Date:
Name of procedure (issue or protocol):
Ask PWPL to invite all involved into the intent.
Pertinent info or details:
Breadth of Intent (1-10) and/or states (1-10):
Procedure type (1 or more): Group I/D - Family I/D - Family Line I/D - Culture I/D
Tools needed:
Prep steps:, support(s):, cleansing(s):, block(s) to be healed:
Building blocks to be installed:
After all steps taken,% of PWPL present.
If 100%, help client put the whole procedure into their circuit with PWPL.
Target work steps:
After all steps are taken,% of PWPL present.
If 100%, help client put the whole procedure into their circuit with PWPL.
(I use a Box metaphor to contain all the procedure remains to give to PWPL.)
If PWPL is not 100%, are there Soul Blocks ? How many?
Should Soul Building Components (SBC s) be created and installed? How many?
Release all Soul Blocks and/or add new Soul Building Components with appropriate level of
intent.
After all Soul Blocks are released and/or SBC s added,% of PWPL present.
If 100%, help client get the whole procedure into their circuit with PWPL.
If PWPL is not 100%, are there Weaknesses ? How many?
Resolve all Weaknesses with the appropriate level of intent.
After all Weaknesses resolved,% of PWPL present.
If 100%, help client put the whole procedure into their circuit with PWPL.
If PWPL is not 100%, are Quality levels 100%? How many levels need to be raised?
Should Latent Qualities for a Procedure type be activated or installed?
Raise all Quality levels and/or activate/install Latent Qualities with appropriate level of intent.
reaso an Quality levels and/or decivate/instant Laterit Qualities with appropriate level of intente.
After all Quality and Latent Quality levels are 100%,% of PWPL present.
If 100%, help client put the whole procedure into their circuit with PWPL.
If PWPL is not 100%, do Support or Ancestral Support circuits need to be made?
How many?
Make all Support and/or Ancestral Support circuits with the appropriate level of intent.
When PWPL is 100%, Quality levels are present, active and 100%, all Support and Ancestral Support circuits are made, help the client put the whole procedure in their circuit
with PWPL.
Check how many more procedures needed for the purpose of work and when: